

Production protocol Baclofen suspension (2-10 mg/mL) with Liquifine® pH4 liquid 500 mL



Preparation instructions

1. Weigh the required quantity of baclofen.
2. Reduce the particles to a fine powder by the use of a mortar and vessel, if necessary.
3. Humidify the powder with a limited quantity of Liquifine® pH4 liquid and homogenize.
4. Add Liquifine® pH4 liquid following the geometric dilution method, to the desired volume.

Production report

Preparation date	/ /	Galenic type	Suspension	
Batch number		Volume		
Raw material	Number analysis register	Required Quantity		Weighed quantity
Baclofen		g		g
Liquifine® pH4 liquid		Ad	mL	Ad mL

Tarra

- _____

Quality control

- Homogeneity: _____

Stability

- 180 days

Storage and labelling

- Protected from light – Room temperature
- Shake before use

Name and Signature Compounder	Name and Signature Pharmacist

MAGIS PHARMA

Keizershoek 336, 2550 Kontich, T: +32 3 457 11 76, F: +32 3 457 43 40
info@magis-pharma.be, www.magis-pharma.be