

# Production protocol

## Hydrocortisone suspension

(1 mg/mL) with Liquifine® pH4 liquid 500 mL



### Preparation instructions

1. Weigh the required quantity of hydrocortisone.
2. Reduce the particles to a fine powder by the use of a mortar and vessel, if necessary.
3. Humidify the powder with a limited quantity of Liquifine® pH4 liquid and homogenize.
4. Add Liquifine® pH4 liquid, following the geometric dilution method, to the desired volume.

### Production report

Preparation date	/ /	Galenic type	Suspension	
Batch number		Volume		
Raw material	Number analysis register	Required Quantity		Weighed quantity
Hydrocortisonum		g		g
Liquifine® pH4 liquid		Ad	mL	Ad mL

### Tarra

- \_\_\_\_\_

### Quality control

- Homogeneity: \_\_\_\_\_

### Stability

- 180 days

### Storage and labelling

- Protected from light – Room temperature
- Shake before use

Name and Signature Compounder	Name and Signature Pharmacist

### MAGIS PHARMA

Keizershoek 336, 2550 Kontich, T: +32 3 457 11 76, F: +32 3 457 43 40  
info@magis-pharma.be, www.magis-pharma.be